

# Multi Para Patient Monitor ECO

Five standard parameters:  
ECG, RESP, NIBP, SPO2, 1-TEMP

Premium Built Quality



12.1" high-resolution color  
TFT display

*“Comprehensive and Powerful  
Monitoring Solutions”*

## Specification:

- SpO2 pitch tone variation
- Multi-display: standard, large font, graphical & Telular trend
- 7-lead ECG waveforms display simultaneously
- Large tabular and graphic trends information storage
- Built-in rechargeable lithium battery
- Supporting for euthymia analysis & drug dose calculation
- Optional parameter Touch Screen, Recorder, 2-IBP, Nellcor, CMS, Trolley, SPO2, Wall mount, VSB port

## ACER BIOMEDICALS LLP

Office: Industrial Unit-80, Block-G, Sector-63, Noida, UP - 201301  
(India) Ph.: +91-120-4166921, E-mail: [info@acerbiomed.com](mailto:info@acerbiomed.com)  
Web.: [www.acerbiomedicals.com](http://www.acerbiomedicals.com), [www.acerbiomed.com](http://www.acerbiomed.com)

**Service Helpline**  
**18008894232**

# Vital Signs Show Perfectly



Trend Graphs



Big Font Interface



OXI-CRG Display



Durable Probes Connector with high quality



Max 8 waveform display simultaneously, max 6 leads ECG waveform display

## Alarming-More safe and reliable

- 3-level audible and visual alarming, Dual alarm light for physiological and technical alarming.

## Data Storage-Stable and Fast

- 400 groups NIBP measurement
- 72 hours Trend data and Trend graphs
- 200 groups alarm events reviewing
- 12 hours wave form reviewing
- Supporting for Arrhythmia analysis and drug dose calculation.

## Product Feature

- Flexible, Light , Easy Operating
- Colourful and Clear 12.1" TFT screen, Blacklight Buttons Available
- Various Language Operating systems
- Built-in detachable rechargeable lithium battery, 2.5 hours battery life